

Arm Positioning – Innovative Solutions



ArmGuard[®] Patient Arm Protector

- **Secure arms without tucking**
Easy access to ports and lines
- **Pad bony prominences**
Foam pads with adjustable straps protect patient, allow immediate access to full arm in an emergency, and put you in control of arm compression
- **Accommodate bariatric patients**
Width-adjustable to accommodate a broad range of patients; accessory pads extend table mattress width
- **Add stability without interference**
Low-profile for robotic procedures and staff comfort, rail-mounted sleds add support in lateral tilt, up to 15°

Positioning Solutions for Arms



ArmGuard® Classic Sled Kit

Includes four sleds: set of two standard and set of two bariatric sleds, which slide under O.R. table mattress.

Item #40004 includes two Standard and two Bariatric Sleds

Requires:

#40412 ArmGuard STANDARD Packs or

#40101 ArmGuard BASIC Packs



ArmGuard® Rail-Mounted Sled

Rail-mounted, width-adjustable sled attaches to O.R. table rail and adds up to 6 in (15 cm) of width.

Item #40018 includes one width-adjustable Rail-Mounted Sled

Requires:

#40412 ArmGuard STANDARD Packs or

#40101 ArmGuard BASIC Packs

Rail clamp (ask about available options)




ArmGuard® Packs

Single-use pads protect arms, hands, and fingers while providing access to ports and lines. Vacuum-packed for effective storage. For use with all ArmGuard sleds.

Case of 12 packs (each contain one left and one right pad):

#40412 ArmGuard STANDARD Packs – pads 24 in (61 cm) long

#40101 ArmGuard BASIC Packs – pads 20 in (51 cm) long



ArmGuard® Width Extenders (optional)

Single-use bariatric accessory pads for high-BMI patients. Available in two sizes to fit mattress thickness and support patient arms. For use with all ArmGuard sleds.

Case of 12 packs (each contains two extender pads):

#40206 ArmGuard Width Extender – mattresses up to 3 in (7.6 cm)

#40207 ArmGuard XL Width Extender – mattresses up to 4 in (10.2 cm)



www.da-surgical.com

info@da-surgical.com

800.261.9953

CE MADE IN USA Pat. www.da-surgical.com/patents ©2024 D. A. Surgical L11143GA7 Rev5