



# Protect Patient – Allow Easy Airway Access



## FaceGuard<sup>®</sup> Head Shield / Instrument Tray

- **Protects patient face**  
Ensures patient protection while providing a convenient surface for surgical instruments
- **Allows easy airway access**  
Adjustable, sliding instrument tray allows unobstructed view and can be quickly removed in an emergency
- **Fits all tables**  
Width-adjustable frame posts
- **Dual-Angle or Flat tray**  
Flat tray for supine and prone procedures; Dual-Angle tray for Trendelenburg

# Compatible with All O.R. Tables



## FaceGuard® Dual-Angle

Dual-Angle tray is designed to provide a flat working surface in supine or Trendelenburg positions.

Item #60701

Requires two rail clamps (ask about available options)



## FaceGuard® Flat

Flat tray is designed for supine and prone procedures in cardiac, orthopedics, and other specialties.

Item #60702

Requires two rail clamps (ask about available options)



Trays can be removed and frame flipped vertically for full access to patient's face in an emergency.



## FaceGuard® Combo

Includes both Dual-Angle and Flat trays, interchangeable on universal frame. Perfect for multi-specialty O.R. suites.

Item #60703

Requires two rail clamps (ask about available options)