


# Stop Patient Sliding – Reduce Cost & Waste



40° Trendelenburg shown with TrenGuard 450 CLASSIC.  
TrenGuard Reusable Lateral Stabilizing Pillow ghosted for visualization of bolster.

## TrenGuard® Trendelenburg Patient Restraint CLASSIC

- **No sliding**  
Patented bolster keeps patients up to 550 lb (250 kg) from sliding
- **No shoulder braces**  
Soft foam lateral stabilizing pillows control body mass shift
- **No chest straps**  
Promotes normal ventilation
- **Efficient and compatible**  
Accommodates under-patient devices and simplifies patient repositioning



CE MADE IN USA

[www.da-surgical.com](http://www.da-surgical.com)

[info@da-surgical.com](mailto:info@da-surgical.com)


800.261.9953

\*compared to viscoelastic static/friction foam pad, shape conforming foam, gel pad, egg crate foam & surgical tape

Pat. [www.da-surgical.com/patents](http://www.da-surgical.com/patents)

©2024 D. A. Surgical L11146GA7 Rev6

# Single-Use Components for Convenience




## TrenGuard® Patient Support Frame

Secures to O.R. table rails with wide-based clamps to support patients up to 550 lb (250 kg) in extreme Trendelenburg up to 40°.


Item #56100

Requires two rail clamps (ask about available options)

Also available: Bariatric and Dynamic frames (ask for details).



450 Wedge



600 Wedge

## TrenGuard® CLASSIC Wedge Packs

Ideal for high-BMI patients and those without discernible cervical curve. Each pack contains wedge bolster (see above), head stabilizing pillow, and lateral stabilizing pillows.

Case of 12:

#55301 TrenGuard 450 CLASSIC Wedge Packs

#56302 TrenGuard 600 CLASSIC Wedge Packs

## TrenGuard® Solo Wedge Packs

Helpful for last-minute substitution. Each pack contains one wedge bolster.

Case of 12:

#54612 TrenGuard 450 Solo Wedge Packs

#54812 TrenGuard 600 Solo Wedge Packs



## TrenGuard® CLASSIC Packs

Disposable components attach to TrenGuard Patient Support Frame with maximum-strength hook and loop fastener. Each pack is vacuum-packed for effective storage and contains:

- A** Patented bolster – prevents patient sliding
- B** Soft, head stabilizing pillow – controls side-to-side rocking of the head
- C** Soft, lateral stabilizing pillows – control body mass shift during transition to Trendelenburg

Available in two sizes:

- 1) **TrenGuard 450** for patients up to 450 lb (204 kg)
- 2) **TrenGuard 600** for patients up to 550 lb (250 kg)

Case of 12:

#55101 TrenGuard 450 CLASSIC Packs

#56101 TrenGuard 600 CLASSIC Packs

Also available: **HYBRID** option, with reusable lateral stabilizing pillows for use with 2-piece single-use packs, available with standard bolster or wedge.



www.da-surgical.com

info@da-surgical.com

800.261.9953

CE MADE IN USA

Pat. www.da-surgical.com/patents

©2024 D. A. Surgical L11146GA7 Rev6